



Mitotic Counting Reproducibility/Feature Study With Pathologists

Brandon D. Gallas PhD

FDA/CDRH/OSEL/DIDSR

Charles Halsey, DVM, PhD and R. Mark Simpson, DVM, PhD

Molecular Pathology Unit, Laboratory of Cancer Biology and Genetics, NCI




Outline

- Motivation
- Study Design and Analysis Concepts
 - Show simulated data
- NIH clinical study
 - Data collection complete
- Reproducibility study
 - Proposal, looking for collaborators and sites

Motivation


- Evaluate new imaging technology
 - Truth not known
 - Reference for comparison determined by readers
 - Intra-reader agreement
 - Inter-reader agreement
 - ***“Reproducibility” Study***
 - Example tasks
 - Impact of viewing mode on mitotic count (digital vs. optical)
 - Impact of stain on mitotic count (H&E vs. PHH3)
- Study designs and analysis methods are lacking
 - Want to account for reader and case variability

Simulated Data




- "Counts" vs. "Counts"
- Scatter plot **shows** agreement

“Correlation” = 0.804744



Simulated Data

- “Counts” vs. “Counts”
- Scatter plot **shows** agreement
- Summarize agreement
 - Correlation is one of many agreement measures




Simulated Data

- Multiple readers evaluate Multiple cases
- Variability from cases
- Variability from readers


Simulated Data

- Intra-Reader
- Intra-modality
 - Requires replicated reads



		Reference Modality, A				
		Dr. 1	Dr. 2	Dr. 3	...	Dr. 8
Reference Modality, A	Dr. 1	Q_1^{AA*}	p_{12}^{AA}	p_{13}^{AA}	...	p_{18}^{AA}
	Dr. 2	p_{21}^{AA}	Q_2^{AA*}	p_{23}^{AA}	...	p_{28}^{AA}
	Dr. 3	p_{31}^{AA}	p_{32}^{AA}	Q_3^{AA*}	...	p_{38}^{AA}

	Dr. 8	p_{81}^{AA}	p_{82}^{AA}	p_{83}^{AA}	...	Q_8^{AA*}



Simulated Data

- Inter-Reader
- Intra-modality
 - No replicated reads

Intra-Modality Agreement		Reference Modality, A				
		Dr. 1	Dr. 2	Dr. 3	...	Dr. 8
Reference Modality, A	Dr. 1	Q_1^{AA*}	p_{12}^{AA}	p_{13}^{AA}	...	p_{18}^{AA}
	Dr. 2	p_{21}^{AA}	Q_2^{AA*}	p_{23}^{AA}	...	p_{28}^{AA}
	Dr. 3	p_{31}^{AA}	p_{32}^{AA}	Q_3^{AA*}	...	p_{38}^{AA}

	Dr. 8	p_{81}^{AA}	p_{82}^{AA}	p_{83}^{AA}	...	Q_8^{AA*}

Study Design & Analysis

- Intra-Reader or Inter-Reader
- Inter-modality
 - No replicated reads
- Examples
 - Optical vs. Digital
 - H&E vs. PHH3

Inter-Modality Agreement		Reference Modality, A				
		Dr. 1	Dr. 2	Dr. 3	...	Dr. 8
New Modality, B	Dr. 1	Q_1^{AB}	P_{12}^{AB}	P_{13}^{AB}	...	P_{18}^{AB}
	Dr. 2	P_{21}^{AB}	Q_2^{AB}	P_{23}^{AB}	...	P_{28}^{AB}
	Dr. 3	P_{31}^{AB}	P_{32}^{AB}	Q_3^{AB}	...	P_{38}^{AB}

	Dr. 8	P_{81}^{AB}	P_{82}^{AB}	P_{83}^{AB}	...	Q_8^{AB}

Intra-Modality Agreement		Reference Modality, A				
		Dr. 1	Dr. 2	Dr. 3	...	Dr. 8
Reference Modality, A	Dr. 1	$Q_1^{AA^*}$	P_{12}^{AA}	P_{13}^{AA}	...	P_{18}^{AA}
	Dr. 2	P_{21}^{AA}	$Q_2^{AA^*}$	P_{23}^{AA}	...	P_{28}^{AA}
	Dr. 3	P_{31}^{AA}	P_{32}^{AA}	$Q_3^{AA^*}$...	P_{38}^{AA}

	Dr. 8	P_{81}^{AA}	P_{82}^{AA}	P_{83}^{AA}	...	$Q_8^{AA^*}$

Study Design & Analysis

New Device Performance

- Is inter-modality agreement worse than intra-modality agreement?

Baseline Performance

New Device Performance

- Is intra-modality agreement worse than intra-modality agreement?

Baseline Performance

Inter-Modality Agreement		Reference Modality, A				
		Dr. 1	Dr. 2	Dr. 3	...	Dr. 8
New Modality, B	Dr. 1	Q_1^{AB}	p_{12}^{AB}	p_{13}^{AB}	...	p_{18}^{AB}
	Dr. 2	p_{21}^{AB}	Q_2^{AB}	p_{23}^{AB}	...	p_{28}^{AB}
	Dr. 3	p_{31}^{AB}	p_{32}^{AB}	Q_3^{AB}	...	p_{38}^{AB}

	Dr. 8	p_{81}^{AB}	p_{82}^{AB}	p_{83}^{AB}	...	Q_8^{AB}

Intra-Modality Agreement		Reference Modality, A				
		Dr. 1	Dr. 2	Dr. 3	...	Dr. 8
Reference Modality, A	Dr. 1	$Q_1^{AA^*}$	p_{12}^{AA}	p_{13}^{AA}	...	p_{18}^{AA}
	Dr. 2	p_{21}^{AA}	$Q_2^{AA^*}$	p_{23}^{AA}	...	p_{28}^{AA}
	Dr. 3	p_{31}^{AA}	p_{32}^{AA}	$Q_3^{AA^*}$...	p_{38}^{AA}

	Dr. 8	p_{81}^{AA}	p_{82}^{AA}	p_{83}^{AA}	...	$Q_8^{AA^*}$

requires replicated readings

Study Design & Analysis

- Intra-reader agreement
 - Summarize with average over all readers
- Inter-reader agreement
 - Summarize with average over all pairs of readers

Two options to produce study endpoints to answer study question

Inter-Modality Agreement		Reference Modality, A				
		Dr. 1	Dr. 2	Dr. 3	...	Dr. 8
New Modality, B	Dr. 1	Q_1^{AB}	p_{12}^{AB}	p_{13}^{AB}	...	p_{18}^{AB}
	Dr. 2	p_{21}^{AB}	Q_2^{AB}	p_{23}^{AB}	...	p_{28}^{AB}
	Dr. 3	p_{31}^{AB}	p_{32}^{AB}	Q_3^{AB}	...	p_{38}^{AB}

	Dr. 8	p_{81}^{AB}	p_{82}^{AB}	p_{83}^{AB}	...	Q_8^{AB}

Intra-Modality Agreement		Reference Modality, A				
		Dr. 1	Dr. 2	Dr. 3	...	Dr. 8
Reference Modality, A	Dr. 1	Q_1^{AA}	p_{12}^{AA}	p_{13}^{AA}	...	p_{18}^{AA}
	Dr. 2	p_{21}^{AA}	Q_2^{AA}	p_{23}^{AA}	...	p_{28}^{AA}
	Dr. 3	p_{31}^{AA}	p_{32}^{AA}	Q_3^{AA}	...	p_{38}^{AA}

	Dr. 8	p_{81}^{AA}	p_{82}^{AA}	p_{83}^{AA}	...	Q_8^{AA}

requires replicated readings

Study Design & Analysis

- Intra-reader agreement
 - Summarize with average over all readers
- Inter-reader agreement
 - Summarize with average over all pairs of readers

Two options to produce study endpoints to answer study question.

- Which best?
- Can do both.

- Pros and Cons for each. Still learning practical and statistical efficiencies.
- Intra-reader requires replicated reading.
- Inter-reader agreement more variable, but averaging over more observations (pairs of readers).

Study Design & Analysis

- What about error bars for study endpoints?
- **MRMC** analysis: account for variability from
 - Multiple Readers
 - Multiple Cases
- Gallas, Anam, Chen, Wunderlich, Zhang (2016), "MRMC analysis of agreement studies." *In Proc. SPIE, 9787, pp. 97870F-97870F-12.*
 - *U-statistics*
 - *Agreement by concordance*
 - *Novel, builds on methods for area under the ROC curve*
 - *Peer reviewed manuscript in development*

Happy to discuss later.


NIH Clinical Study

- PI's from NIH
 - Mark Simpson, DVM, PhD
 - Charles Halsey, DVM, PhD
- Collaborators from FDA: Study Design and Analysis
 - Brandon Gallas, PhD
 - Weijie Chen, PhD
 - Zhiwei Zhang, PhD

NIH Clinical Study


Study Design

- Feature and task:
 - Mitotic Counting
- 12 readers:
 - veterinary pathologists at 4 sites
 - NCI, 2 academic centers, one regional reference lab
- 113 cases:
 - 1 case = 1 slide
 - specimens from canines with oral melanoma
 - **H&E vs. Phosphohistone H3 (PHH3)**



NIH Clinical Study


H&E vs. Phosphohistone H3 (PHH3)



NIH Clinical Study

Study Design

- Counts collected according to standard of care
 - Microscope mode: each pathologist selects 10 evaluation FOVs per case and then counts at 40x
 - Digital mode: simulate standard of care
 - FOV locations saved for each pathologist in digital mode



NIH Clinical Study

Aims

- Compare prognostic ability of mitotic count from
 - H&E Optical, 40x
 - H&E Digital, 20x
 - H&E Digital, 40x
 - PHH3 Optical, 40x
 - PHH3 Digital, 20x
 - PHH3 Digital, 40x

6 modalities = 6 viewing modes

Ultimately, we did not study all modalities and comparisons.
- Inter-reader and Intra-reader agreement analyses
 - Inter-modality and Intra-modality

NIH Clinical Study

Study Design

- Split-plot study design
 - Reduce the total number of reads
 - Reduce workload of each reader
 - Each case read by each reader 3 times instead of 6*
 - Washout time was minimum 1 week*
 - Allow for slides to be shipped

Complicated.
Happy to discuss later.

- Obuchowski, N.; Gallas, B. D. & Hillis, S. L. (2012), 'Multi-Reader ROC Studies with Split-Plot Designs: A Comparison of Statistical Methods.' *Acad Radiol*, **19**, (12), 1508-1517.

NIH Clinical Study

Primary Aim

- Compare prognostic ability of mitotic count
 - H&E slides on a microscope
 - PHH3 slides on a microscope
 - Truth is survival data with some censoring
 - Prognostic ability
= concordance between counts and survival
 - Test non-inferiority with possible superiority

NIH Clinical Study

Preliminary Results

- Sorry. Nothing to share yet.
- Data checks still in progress.

Proposed Reproducibility Study

- PI: Brandon Gallas
- NIH collaborators
 - Mark Simpson, DVM, MD
 - Charles Halsey, DVM, MD
- YOU!?

Proposed Reproducibility Study

Synopsis

- Complement NIH clinical study:
 - mitotic counting
 - H&E vs. PHH3
 - Optical vs. Digital
- Inter-reader and Intra-reader agreement analyses
 - Inter-modality and Intra-modality
- **Study Goal:**
Eliminating location variability yields faster and more precise data

Proposed Reproducibility Study

Synopsis

- Difference from NIH clinical study
 - All pathologists evaluate same pre-selected FOVs
 - FOVs are smaller to fit on a screen
 - Equivalent in area to 10 traditional microscope FOVs
 - Use a reticle on the microscope to get same area as screen
 - Enabled by **eeDAP**:
evaluation environment for digital and analog pathology

Proposed Reproducibility Study

Status


- Protocol under development
- NIH clinical study
 - Source for slides and images
 - Inform study design: modalities, intra- vs. inter-reader
 - Inform study size: readers and cases
 - Results provide basis for comparison
- Collaborators/readers/sites being recruited
 - Hope to involve WSI WG members
 - Need feedback on study protocol
 - Loan and install eeDAP system at your site
 - Find a balance between workload and statistical power (more readers = smaller workload)

Proposed Reproducibility Study

eeDAP

- evaluation environment for digital and analog pathology

- Camera
- Microscope
- Moving stage with multiple slides
- Stage controller
- Joystick for stage control
- Computer and monitor not shown



Brandon D. Gallas, PhD, FDA/CDRH/OSEL/DIDSR

Proposed Reproducibility Study

eeDAP: Register WSI to Glass Slide (Stage)

The screenshot displays the 'Slide Registration' software interface. On the left, a vertical panel contains various controls, with eight numbered callouts (1-8) pointing to specific features:

- 1: Take camera images and stage positions (points to 'Take Stage Position 1' button)
- 2: Take WSI positions (points to 'Take WSI Position 1' button)
- 3: Show camera image (points to 'Camera Image' tab in the inset window)
- 4: Abort experiment (points to 'Abort' button)
- 5: Use last time experiment result (points to 'Load last calibration' button)
- 6: Redo registration (points to 'Refine Registration' button)
- 7: Skip high resolution registration (points to 'Finished Registration' button)
- 8: Show the stage and WSI position (points to 'WSI Position 1' section)

The main window shows a large purple histological slide image with a coordinate grid (0-1000 on both axes). A small green box labeled 'camera ROI' is visible in the top left corner of the slide. An inset window titled 'Registration' is open on the right, showing two side-by-side images: 'Camera Image' (orange background with a red histological section) and 'WSI Image' (white background).

At the bottom of the main window, the text reads: Brandon D. Gallas, PhD, FDA/CDRH/OSEL/DIDSR


Proposed Reproducibility Study

eeDAP: Register WSI to Glass Slide (Stage)

The screenshot displays the 'Slide Registration' software interface. On the left, a vertical panel contains controls for taking camera and WSI coordinates, showing images, and refining registration. Eight numbered callouts point to specific features: 1. 'Take camera images and stage positions' points to the 'Take Stage Position 1' button. 2. 'Take WSI positions' points to the 'Take WSI Position 1' button. 3. 'Show camera image' points to the 'Camera' button. 4. 'Abort experiment' points to the 'Abort' button. 5. 'Use last time experiment result' points to the 'Load last calibration' button. 6. 'Redo registration' points to the 'Refine Registration' button. 7. 'Skip high resolution registration' points to the 'Finished Registration' button. 8. 'Show the stage and WSI position' points to the 'WSI Position 1' section. The main window shows a large purple histological image with a green box indicating a 'camera ROI'. A 'Registration' window on the right compares the 'Camera Image' (orange background) and the 'WSI Image' (purple background), showing a close match in the tissue structure. At the bottom, the text 'Brandon D. Gallas, PhD, FDA/CDRH/OSEL/DIDSR' is visible.

Proposed Reproducibility Study

eeDAP: Collect Data in WSI mode or Microscope mode



Camera image
used to verify
registration during
the study on
microscope.

NOT FOR
PATHOLOGIST
EVALUATION

Parting Messages

Talk summary:

- Presented study design and analysis methods
- Generalize to new readers and new cases
 - = Account for reader and case variability
- Outlined NIH study
- Proposed complementary eeDAP study
 - *Looking for collaborators/sites/readers*
 - *Offer to loan and install an eeDAP system*

Parting Messages

- eeDAP available at <https://github.com/DIDSR/eeDAP>
- Reference
 - Gallas, B. D.; Gavrielides, M. A.; Conway, C.; Ivansky, A.; Keay, T.; Cheng, W.-C.; Hipp, J. & Hewitt, S. M. (2014), 'Evaluation Environment for Digital and Analog Pathology (eeDAP): a platform for validation studies.' *J Med Img*, **1**, (3), **037501**.

Parting Messages

Studies of reproducibility can be a bridge/link between ...

- Technical performance
- Clinical performance
- For example: Track human performance on clinically related task as we vary IQ parameter.

Reductionist approach ...

- Pathologist expertise, case complexity, and clinical context are confounders to technology evaluation and system optimization.
- Limit case variability by defining a narrow task
 - Detect, count, characterize a single feature
- Limit case variability by defining a narrow study population